



ZARC

Apical Shaper[®]

**MORE APICAL PREPARATION
WITH MAXIMUM FLEXIBILITY**

COMPLEMENTARY TO ANY
SYSTEM ON THE MARKET



ApicalShaper® is compatible with all file systems on the market

✓ BlueShaper®

✓ excalibur®

✓ SlimShaper®

Also compatible with:

Reciproc Blue* by VDW

Waveone Gold* by Maillefer

Protaper Gold* by Maillefer

etc...



* Reciproc Blue, Waveone Gold and Protaper Gold are not trademarks of Zarc4Endo.



FROM ODONTOLOGISTS
TO ODONTOLOGISTS

GREATER BACTERIAL REDUCTION MAINTAIN

Apical finishing file with only 3% taper for greater preservation of pericervical dentin



MAXIMUM FLEXIBILITY



GREAT CUTTING EFFICIENCY



**FOUR APICAL PREPARATION OPTIONS:
30, 35, 40 AND 50**



MAXIMUM SECURITY




Apical Shaper®



José Aranguren

ApicalShaper® allows greater cleaning and shaping at the apical level without sacrificing pericervical dentin. We can make very conservative preparations with great cleaning of the canal.



Stéphane Simon

ApicalShaper® is a very useful, simple to use file for shaping large foramen. Its flexibility is incredible.



Valter Fernandes

What I like most of this system is the possibility of being conservative in pericervical dentin without compromising apical disinfection.



MAINTAINING THE ORIGINAL SHAPE OF THE CANAL

Apical Shaper®

The ApicalShaper® family has 4 instruments so that the specialist can choose the one that best suits their clinical needs:


30/03, 35/03, 40/03 and 50/03

	MFD	TAPER	Ø
30	0,78	3%	0,30
35	0,83	3%	0,35
40	0,88	3%	0,40
50	0,98	3%	0,50

Apical Shaper®

Technical characteristics

400-500 RPM, TORQUE 4



- The Blue alloy provides **maximum flexibility** and **resistance to cyclic fatigue**.
- The centered parallelogram section provides a high **cutting capacity** and **great flexibility**.

System Reference Codes

ApicalShaper® 21 mm

ZAS2103003 Z30 21 mm (6 units) sterile
ZAS2103503 Z35 21 mm (6 units) sterile
ZAS2104003 Z40 21 mm (6 units) sterile
ZAS2105003 Z50 21 mm (6 units) sterile
ZAS21ASST1 Z30-Z50 21 mm (4 units) sterile

ApicalShaper® 25 mm

ZAS2503003 Z30 25 mm (6 units) sterile
ZAS2503503 Z35 25 mm (6 units) sterile
ZAS2504003 Z40 25 mm (6 units) sterile
ZAS2505003 Z50 25 mm (6 units) sterile
ZAS25ASST1 Z30-Z50 25 mm (4 units) sterile

ApicalShaper® 31 mm

ZAS3103003 Z30 31 mm (6 units) sterile
ZAS3103503 Z35 31 mm (6 units) sterile
ZAS3104003 Z40 31 mm (6 units) sterile
ZAS3105003 Z50 31 mm (6 units) sterile
ZAS31ASST1 Z30-Z50 31 mm (4 units) sterile



FROM ODONTOLOGISTS
TO ODONTOLOGISTS.

hola@zarc4endo.com
www.zarc4endo.com



ZARC