

ZARC

BlueShaper PRO[®]

Two heat treatments in one file



DualWire[®]


Exceptional flexibility and cut

BlueShaper PRO® is the first system in the world that combines two heat treatments in a single file, thanks to **DualWire® technology**.



José Aranguren

A much more versatile file, cutting, penetration, and flexibility in a single instrument.




Arman Gazi

Combining flexibility and cutting efficiency is an incredible necessity in our canal preparation. My treatments are now safer and easier thanks to DualWire®.



Carlos Franco

With BlueShaper PRO® files, progression in the canal becomes faster and there is a greater sense of security regarding tip resistance.



The **DualWire®** technology, developed by Zarc, allows **greater versatility** of the instrument, offering a **cutting tip** while maintaining the **same flexibility** as the classic BlueShaper®.

BlueShaper PRO®

The efficiency of the new alloy

Thanks to the **dual heat treatment** of the **BlueShaper PRO®**, it allows a more precise and efficient canal preparation.

It features an enhanced cutting tip thanks to its golden 3 mm.



DualWire®

DualWire® technology

The first file system that achieves the incorporation of **2 heat treatments in a single file**, patented by Zarc4Endo.

Designed for cutting-edge endodontics

BlueShaper PRO[®] has been specifically designed to meet the needs of endodontic professionals, providing a **state-of-the-art tool**.

Exceptional durability and performance

It has been manufactured with **high-quality materials** and subjected to **rigorous manufacturing standards**.

It offers outstanding performance thanks to its **increased canal penetration capacity**, increased **resistance to torsional stress**, and improved **cutting efficiency** while maintaining the high flexibility of BlueShaper[®].

BlueShaper PRO[®]

Elevates the quality of endodontic treatment



Enhanced cutting tip



Great torsional resistance



Reduced cyclic fatigue



Maximized flexibility

	TAPER	Ø
Z1	10-2%	0,14
Z2	10-2%	0,17
Z3	5-6%	0,19
Z4	5-6%	0,25

BlueShaper PRO[®]

The BlueShaper PRO[®] family consists of 8 instruments for the specialist to choose the sequence that best suits each clinical case:

Basic sequence:

Z1, Z2, Z3 y Z4

Complete sequence:

ZX, Z1, Z2, Z3, Z4, Z5, Z6 y Z7

Also available separately:

ZX, Z1, Z2, Z3, Z4, Z5, Z6 y Z7

BlueShaper PRO®

Easy incorporation into the clinical routine

BlueShaper PRO® integrates seamlessly into the clinical routine of dentists, providing **easy and consistent results**.

- Exclusive **DualWire®** technology.
- Unique system with two heat treatments in a single file.
- **Greater versatility**, increased cutting efficiency, torsional stress resistance and penetration of the ducts.
- Cyclic fatigue **resistance**.

Zarc, leaders in innovation

Packaging information

BlueShaper PRO® ZX 19 mm

ZBS19PROZX ZX 19 mm (6 units) sterile

BlueShaper PRO® 21 mm

ZBS21ZPRO1 Z1 21 mm (6 units) sterile
ZBS21ZPRO2 Z2 21 mm (6 units) sterile
ZBS21ZPRO3 Z3 21 mm (6 units) sterile
ZBS21ZPRO4 Z4 21 mm (6 units) sterile
ZBS21ZPRO5 Z5 21 mm (6 units) sterile
ZBS21ZPRO6 Z6 21 mm (6 units) sterile
ZBS21ZPRO7 Z7 21 mm (6 units) sterile
ZBS21APRO6 ZX-Z5 21 mm (6 units) sterile
ZBS21APRO4 Z1-Z4 21 mm (4 units) sterile

BlueShaper PRO® 25 mm

ZBS25ZPRO1 Z1 25 mm (6 units) sterile
ZBS25ZPRO2 Z2 25 mm (6 units) sterile
ZBS25ZPRO3 Z3 25 mm (6 units) sterile
ZBS25ZPRO4 Z4 25 mm (6 units) sterile
ZBS25ZPRO5 Z5 25 mm (6 units) sterile
ZBS25ZPRO6 Z6 25 mm (6 units) sterile
ZBS25ZPRO7 Z7 25 mm (6 units) sterile
ZBS25APRO6 ZX-Z5 25 mm (6 units) sterile
ZBS25APRO4 Z1-Z4 25 mm (4 units) sterile

BlueShaper PRO® 31 mm

ZBS31ZPRO1 Z1 31 mm (6 units) sterile
ZBS31ZPRO2 Z2 31 mm (6 units) sterile
ZBS31ZPRO3 Z3 31 mm (6 units) sterile
ZBS31ZPRO4 Z4 31 mm (6 units) sterile
ZBS31ZPRO5 Z5 31 mm (6 units) sterile
ZBS31ZPRO6 Z6 31 mm (6 units) sterile
ZBS31ZPRO7 Z7 31 mm (6 units) sterile
ZBS31APRO6 ZX-Z5 31 mm (6 units) sterile
ZBS31APRO4 Z1-Z4 31 mm (4 units) sterile

BlueShaper Pro® gutta-percha and paper tips are an integral part of the BlueShaper® line.

Gutta-percha

GP-144 Gutta-percha Z3 (60 units)
GP-145 Gutta-percha Z4 (60 units)
GP-146 Gutta-percha Z5 (60 units)
GP-147 Gutta-percha Z6 (60 units)
GP-148 Gutta-percha Z7 (60 units)

Paper points

PP-144 Paper points Z3 (100 units)
PP-145 Paper points Z4 (100 units)
PP-146 Paper points Z5 (100 units)
PP-147 Paper points Z6 (100 units)
PP-148 Paper points Z7 (100 units)

hola@zarc4endo.com
www.zarc4endo.com



ZARC