

BlueShaper PRO[®]

ZARC

	TAPER	Ø
Z5	6-5%	0,30
Z6	5-4%	0,40
Z7	5-4%	0,50

BlueShaper PRO[®]

Z1-Z4 Basic sequence
(85% of cases)

Z5-Z7 If necessary
(15% of cases)



KFile
#10

Glide Path



Z1 **Z2**
RPM: 500
Torque: 4

Shaping



Z3 **Z4** **Z5** **Z6** **Z7**
RPM: 500
Torque: 4

Finishing



*Irrigate, recapitulate with
a 10 file and re-irrigate
after removing each
rotary file*